

MedInformatix System Requirements

A MedInformatix installation requires a workstation for each user who will access the system and a central server to store and process the data. A large installation may require more than one server. The MedInformatix database will be managed by Microsoft SQL Server, which is a software package that will run solely on the server machine. At least one backup storage device is needed to protect the data in case of system failure such as a network storage device, file share or tape drive. Printers also need to be available to print reports, forms and documents. The following topics will give a more detailed explanation of the equipment and software requirements.

The hardware configurations listed below can be used as a guide to server and workstation requirements for MedInformatix Systems. This should be used as a guide only. The final configuration will depend upon the actual use of the system including tolerated down time in the event of hard drive or server failure.

REMOTE ACCESS

The Customer will provide the Acentec Implementation Project and Support Team appropriate remote access to the system for installation, updates and customer support. Remote access using LogMeIn is required. If access is not available or configured using the specified connection types, this will result in a delay in the implementation timeline and support responsiveness.

SERVERS

Each Server is unique to its site. The specifications will vary greatly depending upon the environment. For instance, a site having 50 users may be able to function with a single server whereas a site with 1000 users will require more. If you have questions regarding a server that you already have, call Acentec to see if it will be sufficient.

An installation requires at least one server computer that will store the centralized data and program files for the network. A server machine will provide two functions for the MedInformatix software. It will store the system files and thus act as a file server. It will also run a database management application that will manage the MedInformatix SQL database. A server that runs SQL database software is called a "database server". On a small site a single server machine will handle both server functions. In larger sites, separate server machines will handle each function. Acentec currently recommends Dell servers for the technical support service, value and reliability.

Notes:

- For a fail-safe environment in the event of a hard disk and or software failure, a Raid 5 installation (for hard disk failures) is recommended.
- Other methods of protection such as mirrored drives, SQL server replication, etc. can be installed. The user should review all options and implement accordingly.
- For Database hard drive configuration assume 1 GB per Physician per year for systems using EMR and Practice Management
- Please review the discussion at the end of this document regarding SQL Server capacity and performance.

5- 25 Users – Server Configuration

Single Server for Network Operating/Apps and SQL Server Database

- Windows Essential Business Server®2008 64Bit (Premium Edition)
- O/S – Microsoft Windows Server 2008 64Bit (Standard Edition)
- SQL Database – Microsoft SQL Server 2008 64Bit (Standard Edition)
- Processor – Quad Core Intel® Xeon® X3363, 2.83GHz, 2x6M Cache, 1333MHz FSB
- Memory – 8 GB RAM or more
- Drives –4 Hard Drives of 500GB (raid 5) with global hot spare.
- Optical Drive: 1 DVD drive
- Backup Device / Network Drive Backup / Others

TERMINAL SERVER

Microsoft suggests using remote terminal servers for front end and backup server topology, which means based on the number of remote users accessing your applications, you need calculate the number of terminal servers that you need to load balance your traffic.

Qty 1 – Terminal Server

- O/S – Microsoft Windows Server 2008 64Bit (Standard Edition)
- Processor – Dual Intel X5550, 2.66Ghz, 8M Cache, Turbo, HT, 1333MHz
- Memory – 8 GB RAM
- Drives – 2 Hard Drives of 73 GB each (RAID 1)

Note: Hard disk spaces at varies at this point by number of users using the system.

FAX SERVER

ActiveFax

- O/S – Microsoft Windows Server 2008
- Processor – Intel X5520, 2.26Ghz, 8M Cache, Turbo, HT, 1066MHz
- Memory – 4 GB of RAM
- Drives – 1 Hard Drive of 146 GB
- BrookTrout 4 – 8 Port Modem Card

Note: This component may also be tied to Application Server

NETWORK

Microsoft and Acentec recommend maximizing network hardware to optimize server and workstation performance. Therein, gigabit LAN cards, routers, and switches provide seamless network and SQL application connectivity. In addition, utilizing wireless-G technologies for wireless network components provides ample connectivity performance for Tablet PCs, handheld devices, and all other wireless tools.

WORKSTATIONS

Each user that will access the MedInformatix system will need a workstation. While workstation configurations will vary according to each user's needs and duties, the following can be used as a general guideline for end-user utilization of the MedInformatix system.

Processor

- Single Processor – Pentium IV (Minimum)
- Single Processor – Pentium Core 2Duo (Recommended)

Memory

- 512 MB RAM (Minimum)
- 2 GB RAM (Recommended)

Drives

- 1 Hard Drive – 80 GB (Minimum)
- 1 Hard Drive – 160 GB (Recommended)
- 1 CD Drive (Minimum)
- 1 CDRW/DVD Drive (Recommended)
- Floppy Drive (Optional)

Network Card

- LAN 10/100 Mbps (Minimum)
- Gigabit LAN 10/100/1000 Mbps (Recommended)

O/S & Software

Microsoft XP Professional SP3 (Minimum)
Windows 7 Professional or higher
Internet Explorer 6.0 or higher
Adobe Acrobat Reader

TABLETS (Stylistic/Tablet PC)

Processor

Single Processor—Intel Core Solo (Minimum)
Single Processor—Intel Core2Duo (Recommended)

Memory

512MB RAM (Minimum)
2GB RAM (Recommended)

Drives

1 Hard Drive – 80 GB (Minimum)
1 Hard Drive – 160 GB (Recommended)
1 CD Drive

Network Card

LAN 10/100 Mbps (Minimum)
Gigabit LAN 10/100/1000 Mbps (Recommended)

O/S & Software

Windows XP Professional, Tablet PC 2005 (Minimum)
Windows 7 Professional for Tablet or Higher
Internet Explorer 6.0 or higher
Adobe Acrobat Reader
Microsoft Office One Note

If there are any questions regarding the particular configuration for workstations for your site, please contact Acentec.

Hardware and Accessories

Tablets

<http://solutions.us.fujitsu.com/www/solutions.shtml?solutions/healthcare>

Scanners

<http://www.fujitsu.com/us/services/computing/peripherals/scanners/medical.html>

Label Makers

http://global.dymo.com/enUS/Segments/Office_Solutions.html

Signature Devices

<http://www.wacom.com/industries/medical.php>